

Product datasheet for **TP724922**

PLEKHA7 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human PLEKHA7 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	S287-P381
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human PLEKHA7 Protein is produced by E.coli expression system and the target gene encoding S287-P381 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	144100
UniProt ID:	Q6IQ23
Synonyms:	PH domain-containing family A member 7;Pleckstrin homology domain-containing family A member 7



[View online »](#)