

Product datasheet for **TP724925**

ARHGAP30 Human Recombinant Protein

Product data:

| | |
|---------------------------------------|--|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant Human Rho GTPase-activating protein 30 |
| Species: | Human |
| Expression cDNA Clone or AA Sequence: | 962-1065 |
| Buffer: | Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4 |
| Note: | MADE TO ORDER Recombinant Human Rho GTPase-activating protein 30 is produced by E.coli expression system and the target gene encoding 962-1065 is expressed. |
| Storage: | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles. |
| Stability: | 12 months from date of despatch |
| Locus ID: | 257106 |
| UniProt ID: | Q7Z6I6 |
| Synonyms: | Rho-type GTPase-activating protein 30; ARHGAP30 |



[View online »](#)