

Product datasheet for **TP724930**

CEP57 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Centrosomal protein of 57 kDa
Species:	Human
Expression cDNA Clone or AA Sequence:	S118-R226
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Centrosomal protein of 57 kDa is produced by E.coli expression system and the target gene encoding S118-R226 is expressed.
Stability:	12 months from date of despatch
Locus ID:	9702
UniProt ID:	Q86XR8



[View online »](#)