

Product datasheet for **TP724934**

ASC2 (PYDC1) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Pyrin domain-containing protein 1
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-A89
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Pyrin domain-containing protein 1 is produced by E.coli expression system and the target gene encoding M1-A89 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	260434
UniProt ID:	Q8WXC3
Synonyms:	PAAD-only protein 1; Pyrin-only protein 1; cellular POP1; cPOP1; PYDC1;ASC2; ASC1; POP1; PYC1
Protein Pathways:	NOD-like receptor signaling pathway



[View online »](#)