

Product datasheet for **TP724937**

PPP1R14B Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Protein phosphatase 1 regulatory subunit 14B
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-K147
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Protein phosphatase 1 regulatory subunit 14B is produced by E.coli expression system and the target gene encoding M1-K147 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	26472
UniProt ID:	Q96C90
Synonyms:	Phospholipase C-beta-3 neighbouring gene protein; PPP1R14B; PLCB3N; PNG



[View online »](#)