

Product datasheet for **TP724938**

SAT2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Thialysine N-epsilon-acetyltransferase
Species:	Human
Expression cDNA Clone or AA Sequence:	1-170
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Thialysine N-epsilon-acetyltransferase is produced by E.coli expression system and the target gene encoding 1-170 is expressed.
Stability:	12 months from date of despatch
Locus ID:	112483
UniProt ID:	Q96F10



[View online »](#)