

Product datasheet for **TP724940**

PFDN5 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Prefoldin subunit 5
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-A154
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Prefoldin subunit 5 is produced by E.coli expression system and the target gene encoding M1-A154 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	5204
UniProt ID:	Q99471
Synonyms:	Myc modulator 1; c-Myc-binding protein Mm-1; PFDN5; MM1; PFD5
Protein Families:	Transcription Factors



[View online »](#)