

Product datasheet for **TP724943**

RBM24 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human RNA-binding protein 24
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-G191
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human RNA-binding protein 24 is produced by E.coli expression system and the target gene encoding M1-G191 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	221662
UniProt ID:	Q9BX46
Synonyms:	RNA-binding motif protein 24; RNA-binding region-containing protein 6; RBM24; RNPC6
Protein Families:	Druggable Genome



[View online »](#)