

Product datasheet for **TP724947**

A1CF Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human APOBEC1 complementation factor
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-R122
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human APOBEC1 complementation factor is produced by E.coli expression system and the target gene encoding M1-R122 is expressed.
Stability:	12 months from date of despatch
Locus ID:	29974
UniProt ID:	Q9NQ94



[View online »](#)