

Product datasheet for **TP724947**

A1CF Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human APOBEC1 complementation factor
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-R122
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human APOBEC1 complementation factor is produced by E.coli expression system and the target gene encoding M1-R122 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	29974
UniProt ID:	Q9NQ94
Synonyms:	APOBEC1-stimulating protein; A1CF; ACF; ASP



[View online »](#)