

Product datasheet for **TP724958**

NFU1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human NFU1 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	G10-P254
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human NFU1 Protein is produced by E.coli expression system and the target gene encoding G10-P254 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	27247
UniProt ID:	Q9UMS0
Synonyms:	NFU1 iron-sulfur cluster scaffold homolog, mitochondrial; HIRA-interacting protein 5; HIRIP5



[View online »](#)