

Product datasheet for **TP724963**

CSDC2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Cold shock domain-containing protein C2
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-S153
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Cold shock domain-containing protein C2 is produced by E.coli expression system and the target gene encoding M1-S153 is expressed.
Stability:	12 months from date of despatch
Locus ID:	27254
UniProt ID:	Q9Y534



[View online »](#)