

Product datasheet for **TP724964**

DGKB Human Recombinant Protein

Product data:

| | |
|---------------------------------------|---|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant Human Diacylglycerol kinase beta/DGKB |
| Species: | Human |
| Expression cDNA Clone or AA Sequence: | K30-C140 |
| Buffer: | Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4 |
| Note: | MADE TO ORDER Recombinant Human Diacylglycerol kinase beta/DGKB is produced by E.coli expression system and the target gene encoding K30-C140 is expressed. |
| Stability: | 12 months from date of despatch |
| Locus ID: | 1607 |
| UniProt ID: | Q9Y6T7 |



[View online »](#)