

Product datasheet for **TP724965**

SCIN Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Scinderin
Species:	Human
Expression cDNA Clone or AA Sequence:	1-357
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Scinderin is produced by E.coli expression system and the target gene encoding 1-357 is expressed.
Stability:	12 months from date of despatch
Locus ID:	85477
UniProt ID:	Q9Y6U3



[View online »](#)