

## Product datasheet for **TP724969**

### **BIN1 Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant Human Myc box-dependent-interacting protein 1
<b>Species:</b>	Human
<b>Expression cDNA Clone or AA Sequence:</b>	P400-P593
<b>Buffer:</b>	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
<b>Note:</b>	MADE TO ORDER   Recombinant Human Myc box-dependent-interacting protein 1 is produced by E.coli expression system and the target gene encoding P400-P593 is expressed.
<b>Stability:</b>	12 months from date of despatch
<b>Locus ID:</b>	274
<b>UniProt ID:</b>	<u><a href="#">O00499</a></u>



[View online »](#)