

## Product datasheet for **TP724970**

### SPTBN2 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Spectrin beta chain, non-erythrocytic 2
Species:	Human
Expression cDNA Clone or AA Sequence:	Q2150-K2390
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Spectrin beta chain, non-erythrocytic 2 is produced by E.coli expression system and the target gene encoding Q2150-K2390 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	6712
UniProt ID:	<a href="#">O15020</a>
Synonyms:	Beta-III spectrin; Spinocerebellar ataxia 5 protein; SPTBN2; SCA5



[View online »](#)