

## Product datasheet for **TP724972**

### Periplakin (PPL) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Periplakin
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-K130
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Periplakin is produced by E.coli expression system and the target gene encoding M1-K130 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	5493
UniProt ID:	<a href="#">O60437</a>
Synonyms:	190 kDa paraneoplastic pemphigus antigen; 195 kDa cornified envelope precursor protein; PPL



[View online »](#)