

Product datasheet for **TP724976**

RTN3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Reticulon-3/RTN3(N-His-SUMO)
Species:	Human
Expression cDNA Clone or AA Sequence:	E750-T863
Tag:	N-His-SUMO
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Reticulon-3/RTN3 is produced by E.coli expression system and the target gene encoding E750-T863 is expressed with a His&SUMO tag at the N-terminus.
Stability:	12 months from date of despatch
Locus ID:	10313
UniProt ID:	O95197



[View online »](#)