

Product datasheet for **TP724978**

CALM1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Calmodulin-2
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-K149
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Calmodulin-2 is produced by E.coli expression system and the target gene encoding M1-K149 is expressed.
Stability:	12 months from date of despatch
Locus ID:	801
UniProt ID:	P62158



[View online »](#)