

Product datasheet for **TP724981**

SARS-CoV-2 Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant SARS-CoV-2 NSP7 (C-His)
Species:	SARS-CoV-2
Expression cDNA Clone or AA Sequence:	Ser1-Gln83
Tag:	C-TEV-8His
Buffer:	Supplied as a 0.2 um filtered solution of 20mM Tris, 150mM NaCl, 10%Glycerol, 4mM MgCl ₂ , pH8.0
Note:	Recombinant SARS-CoV-2 NSP7 is produced by our E.coli expression system and the target gene encoding Ser1-Gln83 is expressed with a 8His tag at the C-terminus.
Stability:	12 months from date of despatch
Summary:	The replication of coronavirus is operated by a set of non-structural proteins (nsps) encoded by open reading frame 1a (ORF1a) and ORF1ab in its genome. NSP7 is a component of the coronavirus replicase polyprotein that forms a replication complex. NSP7 is conserved among all coronavirus. It can form a dimer with NSP8 and subsequently assembles into a large hexadecameric structure, which participates in viral replication by acting as a primase. The NSP12-NSP7-NSP8 subcomplex is found to be the minimal core component for the synthesis of coronavirus RNA. NSP7 is of interest for its potential roles in the transcription and replication of the positive-stranded viral RNA genome.



[View online »](#)