

Product datasheet for **TP724983**

SERPINB4 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Serpin B4
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-P390
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Serpin B4 is produced by E.coli expression system and the target gene encoding M1-P390 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	6318
UniProt ID:	P48594
Synonyms:	Leupin; Peptidase inhibitor 11; PI-11; Squamous cell carcinoma antigen 2; SCCA-2; SERPINB4
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)