

## Product datasheet for **TP724984**

### PC4 (SUB1) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human SUB1/RPO2TC1 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-L127
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human SUB1/RPO2TC1 Protein is produced by E.coli expression system and the target gene encoding M1-L127 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	10923
UniProt ID:	<a href="#">P53999</a>
Synonyms:	Activated RNA polymerase II transcriptional coactivator p15;Positive cofactor 4; PC4; SUB1 homolog; p14
Protein Families:	Transcription Factors



[View online »](#)