

Product datasheet for **TP724985**

ERH Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Enhancer of rudimentary homolog
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-K104
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Enhancer of rudimentary homolog is produced by E.coli expression system and the target gene encoding M1-K104 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	2079
UniProt ID:	P84090
Synonyms:	ERH



[View online »](#)