

Product datasheet for **TP724988**

AKAP12 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human A-kinase anchor protein 12
Species:	Human
Expression cDNA Clone or AA Sequence:	E1531-S1782
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human A-kinase anchor protein 12 is produced by E.coli expression system and the target gene encoding E1531-S1782 is expressed.
Stability:	12 months from date of despatch
Locus ID:	9590
UniProt ID:	Q02952



[View online »](#)