

Product datasheet for **TP724989**

TAP1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Antigen peptide transporter 1/APT1
Species:	Human
Expression cDNA Clone or AA Sequence:	Y525-E808
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Antigen peptide transporter 1/APT1 is produced by E.coli expression system and the target gene encoding Y525-E808 is expressed.
Stability:	12 months from date of despatch
Locus ID:	6890
UniProt ID:	Q03518



[View online »](#)