

## Product datasheet for **TP725002**

### TTC38 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Tetratricopeptide repeat protein 38(N-His-SUMO)
Species:	Human
Expression cDNA Clone or AA Sequence:	P363-M468
Tag:	N-His-SUMO
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Tetratricopeptide repeat protein 38 is produced by E.coli expression system and the target gene encoding P363-M468 is expressed with a His&SUMO tag at the N-terminus.
Stability:	12 months from date of despatch
Locus ID:	55020
UniProt ID:	<a href="#">Q5R3I4</a>



[View online »](#)