

Product datasheet for **TP725022**

ARS2 (SRRT) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Serrate RNA effector molecule homolog
Species:	Human
Expression cDNA Clone or AA Sequence:	L527-F876
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Serrate RNA effector molecule homolog is produced by E.coli expression system and the target gene encoding L527-F876 is expressed.
Stability:	12 months from date of despatch
Locus ID:	51593
UniProt ID:	Q9BXP5



[View online »](#)