

## Product datasheet for **TP725022**

### ARS2 (SRRT) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Serrate RNA effector molecule homolog
Species:	Human
Expression cDNA Clone or AA Sequence:	L527-F876
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human Serrate RNA effector molecule homolog is produced by E.coli expression system and the target gene encoding L527-F876 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	51593
UniProt ID:	<a href="#">Q9BXP5</a>
Synonyms:	Arsenite-resistance protein 2; SRRT; ARS2; ASR2



[View online »](#)