

Product datasheet for **TP725027**

GLIPR2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Golgi-associated PR-1 protein/GAPR-1
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-Y90
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Golgi-associated PR-1 protein/GAPR-1 is produced by E.coli expression system and the target gene encoding M1-Y90 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	152007
UniProt ID:	Q9H4G4
Synonyms:	Golgi-associated plant pathogenesis-related protein 1; Glioma pathogenesis-related protein 2; GliPR 2; GLIPR2



[View online »](#)