

Product datasheet for **TP725030**

EHD2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human EH domain-containing protein 2/EHD2
Species:	Human
Expression cDNA Clone or AA Sequence:	L194-E543
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human EH domain-containing protein 2/EHD2 is produced by E.coli expression system and the target gene encoding L194-E543 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	30846
UniProt ID:	Q9NZN4
Synonyms:	PAST homolog 2
Protein Pathways:	Endocytosis



[View online »](#)