

Product datasheet for **TP725034**

MADM (NRBP1) Human Recombinant Protein

Product data:

| | |
|--|--|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant Human Nuclear receptor-binding protein(N-His-SUMO) |
| Species: | Human |
| Expression cDNA Clone or AA Sequence: | S36-T165 |
| Tag: | N-His-SUMO |
| Buffer: | Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4 |
| Note: | MADE TO ORDER Recombinant Human Nuclear receptor-binding protein is produced by E.coli expression system and the target gene encoding S36-T165 is expressed with a His&SUMO tag at the N-terminus. |
| Stability: | 12 months from date of despatch |
| Locus ID: | 29959 |
| UniProt ID: | Q9UHY1 |



[View online »](#)