

Product datasheet for **TP725043**

FARP1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human FARP1
Species:	Human
Expression cDNA Clone or AA Sequence:	V347-D539
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human FARP1 is produced by E.coli expression system and the target gene encoding V347-D539 is expressed.
Stability:	12 months from date of despatch
Locus ID:	10160
UniProt ID:	Q9Y4F1



[View online »](#)