

Product datasheet for **TP725044**

SRPR beta (SRPRB) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human SR-beta
Species:	Human
Expression cDNA Clone or AA Sequence:	R58-A271
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human SR-beta is produced by E.coli expression system and the target gene encoding R58-A271 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	58477
UniProt ID:	Q9Y5M8
Synonyms:	Signal recognition particle receptor subunit beta; Protein APMCF1; SRPRB
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)