

Product datasheet for **TP725050**

GATD3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human GATD3 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	A42-K268
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human GATD3 Protein is produced by E.coli expression system and the target gene encoding A42-K268 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	8209
UniProt ID:	P0DPI2
Synonyms:	Glutamine amidotransferase-like class 1 domain-containing protein 3, mitochondrial; GATD3A



[View online »](#)