

Product datasheet for **TP725051**

Methylmalonyl Coenzyme A mutase (MUT) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human MCM
Species:	Human
Expression cDNA Clone or AA Sequence:	E451-V750
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human MCM is produced by E.coli expression system and the target gene encoding E451-V750 is expressed.
Stability:	12 months from date of despatch
Locus ID:	4594
UniProt ID:	P22033



[View online »](#)