

Product datasheet for **TP725052**

ECHS1 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Enoyl-CoA hydratase, mitochondrial
Species:	Human
Expression cDNA Clone or AA Sequence:	A28-Q290
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Enoyl-CoA hydratase, mitochondrial is produced by E.coli expression system and the target gene encoding A28-Q290 is expressed.
Stability:	12 months from date of despatch
Locus ID:	1892
UniProt ID:	P30084



[View online »](#)