

## Product datasheet for **TP725053**

### PCBP3 Human Recombinant Protein

#### Product data:

|                                       |  |
|---------------------------------------|--|
| Product Type:                         | Recombinant Proteins   |
| Description:                          | Recombinant Human Poly(rC)-binding protein 3   |
| Species:                              | Human  |
| Expression cDNA Clone or AA Sequence: | M33-L371   |
| Buffer:                               | Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4  |
| Note:                                 | MADE TO ORDER   Recombinant Human Poly(rC)-binding protein 3 is produced by E.coli expression system and the target gene encoding M33-L371 is expressed. |
| Stability:                            | 12 months from date of despatch  |
| Locus ID:                             | 54039  |
| UniProt ID:                           | <a href="#">P57721</a>   |



[View online »](#)