

## Product datasheet for **TP725057**

### SCCPDH Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human SCCPDH Protein(N-His-SUMO)
Species:	Human
Expression cDNA Clone or AA Sequence:	M1-Y274
Tag:	N-His-SUMO
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human SCCPDH Protein is produced by E.coli expression system and the target gene encoding M1-Y274 is expressed with a His&SUMO tag at the N-terminus.
Stability:	12 months from date of despatch
Locus ID:	51097
UniProt ID:	<a href="#">Q8NBX0</a>



[View online »](#)