

Product datasheet for **TP725059**

ECHDC3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human ECHDC3
Species:	Human
Expression cDNA Clone or AA Sequence:	R29-V303
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human ECHDC3 is produced by E.coli expression system and the target gene encoding R29-V303 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	79746
UniProt ID:	Q96DC8
Synonyms:	Enoyl-CoA hydratase domain-containing protein 3, mitochondrial



[View online »](#)