

Product datasheet for **TP725063**

Cytochrome P450 3A4 (CYP3A4) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Cytochrome P450 3A4
Species:	Human
Expression cDNA Clone or AA Sequence:	D154-A503
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Cytochrome P450 3A4 is produced by E.coli expression system and the target gene encoding D154-A503 is expressed.
Stability:	12 months from date of despatch
Locus ID:	1576
UniProt ID:	P08684



[View online »](#)