

Product datasheet for **TP725068**

KCNK10 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Potassium channel subfamily K member 10/KCNK10
Species:	Human
Expression cDNA Clone or AA Sequence:	R379-N538
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Potassium channel subfamily K member 10/KCNK10 is produced by E.coli expression system and the target gene encoding R379-N538 is expressed.
Stability:	12 months from date of despatch
Locus ID:	54207
UniProt ID:	P57789



[View online »](#)