

Product datasheet for **TP725069**

GABRE Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Gamma-aminobutyric acid receptor subunit epsilon(N-His-SUMO)
Species:	Human
Expression cDNA Clone or AA Sequence:	N29-K253
Tag:	N-His-SUMO
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Gamma-aminobutyric acid receptor subunit epsilon is produced by E.coli expression system and the target gene encoding N29-K253 is expressed with a His&SUMO tag at the N-terminus.
Stability:	12 months from date of despatch
Locus ID:	2564
UniProt ID:	P78334



[View online »](#)