

Product datasheet for **TP725070**

CENPE Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Centromere-associated protein E
Species:	Human
Expression cDNA Clone or AA Sequence:	E2401-R2500
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Centromere-associated protein E is produced by E.coli expression system and the target gene encoding E2401-R2500 is expressed.
Stability:	12 months from date of despatch
Locus ID:	1062
UniProt ID:	Q02224



[View online »](#)