

Product datasheet for **TP725071**

KCNH2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human KCNH2 Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	F860-A990
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human KCNH2 Protein is produced by E.coli expression system and the target gene encoding F860-A990 is expressed.
Stability:	12 months from date of despatch
Locus ID:	3757
UniProt ID:	Q12809



[View online »](#)