

Product datasheet for **TP725074**

BRDT Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human Bromodomain testis-specific protein/BRDT
Species:	Human
Expression cDNA Clone or AA Sequence:	Q795-M926
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human Bromodomain testis-specific protein/BRDT is produced by E.coli expression system and the target gene encoding Q795-M926 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	676
UniProt ID:	Q58F21
Synonyms:	Cancer/testis antigen 9; CT9; RING3-like protein
Protein Families:	Protein Kinase, Transcription Factors



[View online »](#)