

Product datasheet for **TP725077**

TRPV3 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human TrpV3
Species:	Human
Expression cDNA Clone or AA Sequence:	N671-V790
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human TrpV3 is produced by E.coli expression system and the target gene encoding N671-V790 is expressed.
Stability:	12 months from date of despatch
Locus ID:	162514
UniProt ID:	Q8NET8



[View online »](#)