

Product datasheet for **TP725079**

TRPC4 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human TrpC4
Species:	Human
Expression cDNA Clone or AA Sequence:	F746-L977
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human TrpC4 is produced by E.coli expression system and the target gene encoding F746-L977 is expressed.
Stability:	12 months from date of despatch
Locus ID:	7223
UniProt ID:	Q9UBN4



[View online »](#)