

## Product datasheet for **TP725080**

### PIK3CD Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human PI3-kinase subunit delta
Species:	Human
Expression cDNA Clone or AA Sequence:	A396-L580
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER   Recombinant Human PI3-kinase subunit delta is produced by E.coli expression system and the target gene encoding A396-L580 is expressed.
Stability:	12 months from date of despatch
Locus ID:	5293
UniProt ID:	<u><a href="#">O00329</a></u>



[View online »](#)