

Product datasheet for **TP725081**

PDE2A Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human PDE2A Protein
Species:	Human
Expression cDNA Clone or AA Sequence:	A18-R226
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human PDE2A Protein is produced by E.coli expression system and the target gene encoding A18-R226 is expressed.
Stability:	12 months from date of despatch
Locus ID:	5138
UniProt ID:	<u>O00408</u>



[View online »](#)