

Product datasheet for **TP725085**

MAP3K9 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant Human MAP3K9
Species:	Human
Expression cDNA Clone or AA Sequence:	R575-P821
Buffer:	Supplied as a 0.2 um filtered solution of 500mM PBS, 300mM NaCl, 10%Glycerol, pH7.4
Note:	MADE TO ORDER Recombinant Human MAP3K9 is produced by E.coli expression system and the target gene encoding R575-P821 is expressed.
Storage:	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening.Please minimize freeze-thaw cycles.
Stability:	12 months from date of despatch
Locus ID:	4293
UniProt ID:	P80192
Synonyms:	Mitogen-activated protein kinase kinase kinase 9; Mixed lineage kinase 1; MLK1; PRKE1
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)